



Goodhead, L., & MacMillan, F. (2017). Measuring osmosis and hemolysis of red blood cells. *Advances in Physiology Education*, 41(2), 298-305. <https://doi.org/10.1152/advan.00083.2016>

Peer reviewed version

License (if available):  
Unspecified

Link to published version (if available):  
[10.1152/advan.00083.2016](https://doi.org/10.1152/advan.00083.2016)

[Link to publication record in Explore Bristol Research](#)  
PDF-document

This is the author accepted manuscript (AAM). The final published version (version of record) is available online via APS at <http://advan.physiology.org/content/41/2/298.long>. Please refer to any applicable terms of use of the publisher.

## University of Bristol - Explore Bristol Research

### General rights

This document is made available in accordance with publisher policies. Please cite only the published version using the reference above. Full terms of use are available:  
<http://www.bristol.ac.uk/red/research-policy/pure/user-guides/ebr-terms/>

1 MEASURING OSMOSIS AND HAEMOLYSIS OF RED BLOOD CELLS

2 Lauren K Goodhead and Frances M MacMillan

3 School of Physiology, Pharmacology and Neuroscience, University of Bristol, University Walk, Bristol, BS8

4 1TD, United Kingdom

5

6 **CORRESPONDING AUTHOR**

7 Lauren K Goodhead

8 School of Physiology, Pharmacology and Neuroscience, University of Bristol, University Walk, Bristol, BS8

9 1TD, United Kingdom

10 Tel: (0117) 3312369

11 Lauren.Goodhead@bristol.ac.uk

12

13 **ABBREVIATED TITLE**

14 Measuring osmosis and haemolysis of red blood cells

15

16 **KEYWORDS**

17 Haematocrit

18 Hanging tissue fluids

19 Osmolarity

20 Tonicity

21

22 **ABSTRACT**

23 Since the discovery of the composition and structure of the mammalian cell membrane, biologists have had  
24 a clearer understanding of how substances enter and exit the cell's interior. The selectively-permeable nature  
25 of the cell membrane allows the movement of some solutes and prevents the movement of others. This has  
26 important consequences for cell volume and the integrity of the cell and, as a result, is of utmost clinical  
27 importance, for example in the administration of isotonic intravenous infusions. The concepts of osmolarity

28 and tonicity are often confused by students as impermeant isosmotic solutes such as NaCl are also isotonic;  
29 however, isosmotic solutes such as urea are actually hypotonic due to the permeant nature of the  
30 membrane. By placing red blood cells in solutions of differing osmolarities and tonicities, this experiment  
31 demonstrates the effects of osmosis and the resultant changes in cell volume. Using haemoglobin standard  
32 solutions, where known concentrations of haemoglobin are produced, the proportion of haemolysis and the  
33 effect of this on resultant haematocrit can be estimated. No change in cell volume occurs in isotonic NaCl,  
34 and by placing blood cells in hypotonic NaCl incomplete haemolysis occurs. By changing the bathing solution  
35 to either distilled water or isosmotic urea, complete haemolysis occurs due to their hypotonic effects. With  
36 the use of animal blood in this practical, students gain useful experience in handling tissue fluids and  
37 calculating dilutions and can appreciate the science behind clinical scenarios.

38

## 39 **INTRODUCTION**

### 40 **Objectives and Overview**

41 The movement of water and small molecules across the selectively-permeable membranes of mammalian  
42 cells is a fundamental concept of physiology. These processes can be difficult for students to visualise and  
43 appreciate and it is often left to images in textbooks or online animations to explain such movements. This  
44 practical uses animal blood bathed in solutions with differing osmolarities and tonicities to explore the  
45 concept of water movement by osmosis and the resultant haemolysis that can occur when red blood cells  
46 are exposed to hypotonic solutions. Students are given the opportunity to handle body fluids, practise  
47 preparing dilutions and make accurate observations.

48

### 49 **Background**

50 In 1925, Gorter and Grendel (6) were the first to report the bilayer nature of the cell membrane. The structure  
51 of the cell membrane was further advanced by the work of Singer and Nicolson (18) who described the  
52 presence and location of proteins in the bilayer and developed the fluid mosaic model. In the mammalian  
53 cell membrane, the phospholipid bilayer alone is permeable to some substances such as oxygen, a small non-  
54 polar molecule, and partially permeable to water but some substances such as charged ions and glucose are

impermeant without the additional presence of protein channels and transporters in the membrane. The combined properties of the phospholipid and proteins has resulted in the use of the term the 'selectively-permeable' membrane (3,9). The extent to which solutes can cross the cell membrane dictates the tonicity of extracellular fluids and therefore the size and shape of cells from the resultant osmotic water movement (19). Knowledge of the structure and function of cell membranes and the movement of substances across the membrane is fundamental to all biomedical science disciplines and is often taught in early parts of undergraduate courses.

Osmosis is the movement of water down its osmotic gradient across a selectively-permeable membrane (5). The establishment of an osmotic pressure gradient, i.e. the pressure required to prevent the movement of water down its gradient, is a result of the difference in numbers of impermeant particles in solution on either side of the membrane (14). Water can move directly through the cell membrane; however, due to the lipid bilayer nature of the membrane this process is relatively slow. It was the discovery of water carrying pore-forming proteins known as aquaporins (16) that helped improve knowledge of how water moves from intracellular to extracellular fluid and vice versa. Water balance is crucial in homeostasis; hormones such as antidiuretic hormone (ADH) and atrial natriuretic peptide (ANP) are released in response to changes in plasma composition and volume respectively, and act on the kidney to regulate plasma osmolarity and volume.

The osmolarity of a solution is determined by the total number of particles present, known as osmolyte particles, and is not affected by the identity of these molecules (19). The higher the osmolarity of a solution, the greater the concentration of osmolytes and the physical properties of a solution such as osmotic pressure and freezing point will be dependent on the concentration of osmolytes in solution. Osmolarity is calculated from the sum of the molar concentration of each solute multiplied by the osmotic coefficient for that solute. The osmotic coefficient is determined by the degree to which a solute (e.g. an ionic compound) dissociates in solution therefore an osmotic coefficient of '1' indicates that the solute completely dissociates in solution.

81 For example, to calculate the osmolarity of a 0.9% weight/volume NaCl (MW 58.44) solution firstly the  
82 molarity is calculated by:

83

84 Molarity of a % w/v solution (M) = % solution in g/litre ÷ molecular mass of the solute

85  $0.154 \text{ M} = 9 \text{ g/litre} \div 58.44 \text{ g/mol}$

86

87 To calculate the osmolarity, given that NaCl dissociates into two ions ( $\text{Na}^+$  and  $\text{Cl}^-$ ) in solution and has an  
88 osmotic coefficient of 0.93, the following equation is used:

89 Osmolarity of solution (osmol/l) = molarity (M) x number of osmoles produced by dissociation x osmotic  
90 coefficient

91  $0.286 \text{ osmol/l or } 286 \text{ mosmol/l} = 0.154 \text{ M} \times 2 \times 0.93$

92

93 Osmolarity and tonicity are often used interchangeably by students but they are not the same. Tonicity refers  
94 to the effect a solution has on cell volume as a result of the permeability of the membrane to that solute.  
95 Tonicity is therefore determined by the osmolarity and whether the solute can cross the cell membrane; it is  
96 the concentration of the impermeant solutes alone that determines tonicity. When comparing fluid  
97 concentrations to that of extracellular body fluid, the terms isotonic, hypertonic and hypotonic are used  
98 rather than osmolarity as they describe the effect the solution has on cell volume which is of physiological  
99 significance. The tonicity will result in: no net movement of water (isotonic), net flow of water out of a cell  
100 (hypertonic), or net flow of water into a cell (hypotonic). Two solutions that are isosmotic may not be isotonic.  
101 A key example is isosmotic urea and isosmotic NaCl. Both urea and NaCl have the same osmolarity, having  
102 the same total number of osmolyte particles, however the membrane is permeable to urea which will freely  
103 diffuse across the cell membrane and impermeable to NaCl. An isosmotic urea is therefore hypotonic  
104 compared to an isosmotic and isotonic solution of the impermeant NaCl. As a result, the volume of a cell is  
105 determined by the solution in which it is being bathed **and** whether the cell's membrane is permeable to the  
106 solute. If a membrane is not equally permeable to all solutes then a difference in water movement will be  
107 observed that is not explained by osmolarity alone and hence an additional term, tonicity, is required.

108 Hypotonic solutions lead to cell swelling and eventual rupture or lysis if the resultant osmotic movement of  
109 water is great enough. In the case of red blood cells this is referred to as haemolysis (4).

110

111 Knowledge of osmosis and tonicity is crucial in understanding the movement of fluids in the body. These  
112 concepts are fundamental in normal physiological processes; one example is that of water reabsorption in  
113 the kidney as increases in osmolarity are detected by the hypothalamus and stimulate the secretion of ADH  
114 resulting in greater water retention and excretion of more concentrated urine (7). Osmosis and tonicity are  
115 important clinically as the failure of the body to respond to changes in osmolarity, or the failure to release  
116 ADH, result in the condition diabetes insipidus. Another important concept is the diagnosis of the different  
117 types of dehydration and the administration of appropriate intravenous fluids (2). In this practical, using easy  
118 to obtain red blood cells as model cells (1), students can explore the concepts of membrane permeability,  
119 osmosis, osmotic pressure, tonicity and haemolysis whilst also learning key laboratory skills such as making  
120 dilution series and handling tissue fluids.

121

## 122 **Learning Objectives**

123 After completing this activity, the student will be able to:

- 124 1. CONTENT KNOWLEDGE: Define key terms used in explaining concentration, osmolarity, osmotic  
125 pressure and tonicity
- 126 2. CONTENT KNOWLEDGE: Calculate the osmolarity of a solution
- 127 3. CONTENT KNOWLEDGE: Describe and explain the consequences of bathing red blood cells in  
128 solutions of differing tonicity
- 129 4. PROCESS SKILLS: Handle mammalian blood samples safely
- 130 5. PROCESS SKILLS: Prepare standard saline solutions
- 131 6. PROCESS SKILLS: Measure haematocrit and estimate haemoglobin concentration
- 132 7. PROCESS SKILLS: Carry out experiments with careful planning, accurate observation and recording of  
133 results

134

135 **Activity Level**

136 This activity is used to teach students in their first year of undergraduate study in physiology. This practical  
137 is used on our Physiological Sciences programme and Veterinary Science programme but would also be  
138 suitable for other biomedical science or healthcare professional programmes such as medicine.

139

140

141 **Prerequisite Student Knowledge or Skills**

142 Before undertaking this activity, students should have a basic understanding of:

143 Homeostasis and the proportions of fluid in intracellular and extracellular compartments

144 The definition of a solute, a solvent and a solution

145 The concept of osmosis and the movement of water across a selectively-permeable membrane

146

147 Students should know how to:

148 Perform basic calculations to work out volumes required for concentrations

149 Use pipettes to create serial dilutions from stock solutions

150 Collect data carefully and accurately

151 Observe safe laboratory practices

152

153 **Time Required**

154 This practical is run in a 3 hour laboratory time-slot. The practical is completed within one session; however,  
155 it is expected that students complete their pre-reading of the laboratory notes which explain the concepts of  
156 osmolarity, tonicity and how to calculate osmolarity (to aid in achieving content learning objectives 1 and 2)  
157 and an online pre-practical quiz before they come to the practical. This preparation work is expected to take  
158 around 1 hour.

159

160 **METHOD**

161 **Equipment and Supplies**

162 The following equipment and supplies are needed:

163 Solutions

164 1. Distilled water (20 ml per pair of students)

165 2. 2.7% w/v NaCl solution (2.7 g NaCl per 100 ml of distilled water) (20 ml per pair of students plus that  
166 required for non-haemolysed blood preparation). This stock solution is used to prepare all other NaCl  
167 solutions in the experiment

168 3. Isosmotic urea solution (17.1 g/l) (5 ml per pair of students plus that required for haemolysed blood  
169 preparation)

170 4. Fresh mammalian blood. This blood is referred to for the rest of the experiment as non-haemolysed  
171 blood. We find that there are no appreciable differences in the outcome of the experiment depending  
172 on which species blood is used although values of haemolysis can vary. Obtaining mammalian blood  
173 supplies can be problematic if obtained locally direct from an abattoir; however, blood can also be  
174 purchased online for example <http://www.rockland-inc.com/blood-products.aspx>. For a class of around  
175 200 students working in pairs, approximately 1.5 L of blood is required (approximately 11 ml blood per  
176 pair of students and allowing extra for repeat experiments if required). The blood must be heparinised  
177 before use to prevent clotting by the addition of heparin sodium (5,000 I.U./ml per 1.5 litres blood). This  
178 blood is then used to produce the haemolysed and non-haemolysed blood as follows:

179 5. Haemolysed blood. To prepare the haemolysed blood in manageable volumes, 250ml of non-  
180 haemolysed blood is measured into a 600ml beaker together with 250 ml urea solution (17.1g/l) and  
181 stirred. The tonicity of the urea and resultant osmotic water movement results in haemolysis of the cells  
182 and this will form the blood used for the production of the haemoglobin standards that will be used to  
183 assess the degree of haemolysis in the experiment. Decant 10 ml of the haemolysed blood into 50  
184 centrifuge tubes (one per pair of students), labelled H for haemolysed blood and centrifuge at 6000 rpm  
185 for 2 minutes. Repeat depending on quantities of blood required i.e. if 1l required repeat once.

186 6. Non-haemolysed blood. To prepare the non-haemolysed blood in manageable volumes, 275 ml of non-  
187 haemolysed blood (from the original heparinised fresh mammalian blood) is prepared by the addition  
188 of 275 ml of 0.9% w/v saline and stirred gently. This forms the non-haemolysed blood which will be used



for the main part of the experiment at an equal concentration to the haemolysed blood. Decant 11 ml of the non-haemolysed blood into 50 centrifuge tubes (one per pair of students) labelled N for non-haemolysed blood. Repeat depending on quantities of blood required i.e. if 1l required repeat once.

An assumption is made that the haemoglobin concentration of the original blood sample is 15 g/dl; but as the haemolysed blood is diluted 1:1 with isosmotic urea (17.1g/l) and the equivalent non-haemolysed blood is diluted 1:1 with isosmotic (0.9% w/v NaCl), the haemoglobin concentration of both blood samples is therefore assumed to be 7.5 g/dl (75 g/litre).

#### Equipment

1. 600 ml glass beakers (2 for blood preparation)
2. 500 ml measuring cylinders (2 for blood preparation)
3. Stirring rods (2 for blood preparation)
4. 25 ml glass beakers for water, 2.7% w/v NaCl and urea distribution (3 per pair of students)
5. 1.5 ml plastic Eppendorf tubes with hinged cap (11 per pair of students)
6. 10 ml plastic centrifuge tubes with cap (10 per pair of students)
7. Centrifuge tube racks (1 per pair of students)
8. 75 µl glass microhaematocrit tubes (Hawksley catalogue no. 01603) (6 per pair of students)
9. Plasticine
10. Centrifuge with centrifuge tube rotor and microhaematocrit tube rotor (Hettich EBA21 centrifuge with 1416 rotor and 1450 haematocrit rotor)
11. Haematocrit readers (Hawksley) or 30 cm rulers (a number of readers/rulers can be shared between pairs of students)
12. 1.5 ml disposable plastic pipettes or equivalent Gilson pipettes if available (3 disposable pipettes per pair of students)
13. Marker pens (1 per pair of students)
14. White paper (1 sheet per pair of students)

#### **Human or Animal Subjects**

216 The animal blood used in this experiment is obtained as a by-product from a local abattoir and therefore the  
217 animals are not slaughtered for the purpose of this experiment.

218

## 219 **Instructions**

### 220 *Preparation prior to the practical*

221 In advance of the class, students must calculate (1) the volume of distilled water and 2.7% w/v NaCl stock  
222 solution required to produce 9 ml each of 0.9 and 0.45% w/v saline solutions; (2) the volumes of haemolysed  
223 blood and 0.9% w/v NaCl (ml) required to produce 1.5 ml % haemoglobin concentrations and (3) the  
224 haemoglobin concentration (g/dl) in results tables 1-3 provided in their lab books.

225

226 In our programmes this is the first practical that students will have had to calculate and make serial dilutions  
227 and handle blood; two key, but challenging, transferable skills.

228

### 229 *Making saline solutions and haemoglobin standards*

230 We recommend that students carry out this practical working in groups of two or three. Students begin the  
231 practical by making a set of standard solutions of haemolysed blood of known haemoglobin concentration  
232 to use later in the experiment which they will compare against the unknown haemoglobin containing  
233 solutions they will produce. Haemolysed blood is used to create these haemoglobin standards as this contains  
234 red blood cells that have already fully lysed in urea and all the haemoglobin has been released into the  
235 solution. The steps below take the students through the practical:

236 1. Using the 2.7% w/v NaCl and 1.5 ml pipettes provided, prepare 9 ml each of 0.9% w/v NaCl and 0.45% w/v  
237 NaCl solutions in two labelled 10ml plastic centrifuge tubes from the dilutions calculated below in Table 1  
238 (values underlined are calculated by the students in advance of the class).

239

### 240 **Table 1: Dilutions calculations for saline solutions**

241

242 2. Using the marker pen, label five Eppendorf tubes – 100%, 66%, 33%, 7% and 1% which will represent the  
243 percentage of haemolysed blood to be added to these Eppendorf tubes. Using the volumes calculated in  
244 Table 2, use 1.5 ml pipettes to add the appropriate volumes of 0.9% w/v NaCl solution and haemolysed blood  
245 to each labelled Eppendorf tube (values underlined are calculated by the students in advance of the class).  
246 Gently invert the tube containing the haemolysed blood before use to ensure the blood is evenly mixed and  
247 once filled also invert each Eppendorf tube to ensure mixing.

248

249 **Table 2: Dilutions calculations for haemoglobin standards from haemolysed blood**

250

251 The calculations performed in Table 3 provide reference haemoglobin concentrations for each haemoglobin  
252 standard (values underlined are calculated by the students in advance of the class).

253

254 **Table 3: Haemoglobin concentration in each standard solution**

255

256 3. Lay out the five mixed haemolysed blood/0.9% w/v NaCl solutions in the Eppendorf tubes on a blank sheet  
257 of white paper to observe the colours. The colours of the haemoglobin standards should range from  
258 translucent pink to translucent red and should look similar to Figure 1. These haemoglobin standards will be  
259 used later in the experiment and should be kept to one side until then.

260

261 **Figure 1: The haemoglobin standards produced from haemolysed blood with 1% haemolysed blood to the**  
262 **left of the image and 100% haemolysed blood on the right. Image courtesy of the University of Bristol.**

263

264 *Investigating the effects of tonicity on red blood cells*

265 The next part of the experiment investigates the effects of membrane-permeable and membrane-  
266 impermeable solutions of differing concentrations on whole red blood cells using the non-haemolysed blood  
267 sample. Both haematocrit and % haemolysis will be estimated. The haematocrit will indicate the degree to  
268 which red blood cells swell or shrink when exposed to the different solutions but does not take into account

269 if haemolysis has occurred. Percentage haemolysis gives a measure of the degree of haemolysis of the  
270 samples and can be used to determine if red blood cells have swollen and burst. Haematocrit alone cannot  
271 distinguish between cell shrinkage and a combination of swelling and lysis.

272 4. Label an additional six 10 ml plastic centrifuge tubes from 1-6. Gently invert the tube containing the non-  
273 haemolysed blood several times before use to ensure an even suspension of red blood cells.

274 5. Prepare the centrifuge tubes as follows:

275 Tube 1 – 1.5 ml non-haemolysed blood + 1.5 ml 2.7% w/v NaCl

276 Tube 2 – 1.5 ml non-haemolysed blood + 1.5 ml 0.9% w/v NaCl

277 Tube 3 – 1.5 ml non-haemolysed blood + 1.5 ml 0.45% w/v NaCl

278 Tube 4 – 1.5 ml non-haemolysed blood + 1.5 ml distilled water

279 Tube 5 – 1.5 ml non-haemolysed blood + 1.5 ml isosmotic urea

280 Tube 6 – 3 ml non-haemolysed blood

281 When filled, gently invert the centrifuge tubes several times to ensure the blood is mixed and leave for 10  
282 minutes before proceeding with the next step.

283 6. The blood solutions are then prepared for centrifuging to allow the measurement of the packed cell volume  
284 (haematocrit) of each sample. Label six glass microhaematocrit tubes 1-6 to correspond to the samples in the  
285 plastic centrifuge tubes. In turn, invert each centrifuge tube several times to ensure even dispersal of red  
286 blood cells and then dip the corresponding microhaematocrit tube in the blood until capillary action has filled  
287 the glass tube. Seal the bottom of the microhaematocrit tube with a small plug of plasticine by twisting the  
288 bottom of the tube in a tray of plasticine.

289 7. Centrifuge the microhaematocrit tubes at 6,000 rpm for 2 minutes using the microhaematocrit tube rotor  
290 until the cells have packed together at the bottom of the tube leaving the fluid (supernatant) above.

291 It is not expected that the centrifuges are operated by the students, in our laboratory students bring their  
292 samples to the shared laboratory centrifuges and these are run by experienced demonstrators or technicians.

293

294 *Measuring haematocrit*

295 8. After centrifuging, measure the haematocrit of each sample using a haematocrit reader and read off the  
296 % haematocrit. If haematocrit readers are difficult to obtain, a ruler can be used instead. By this method,  
297 measure the total length of the column of fluid and the length of the column of packed cells and calculate  
298 the proportion of the total column that is made up of packed cells at the bottom. This percentage is the  
299 haematocrit. Record the haematocrit readings in the *observed haematocrit* column in Table 4.

300

301 **Table 4: Haematocrit measurements for non-haemolysed blood**

302

303 9. With the exception of tube 6, the haematocrit readings measured are for blood diluted 50:50 with a saline  
304 solution. Therefore complete the *dilution factor* column with a dilution factor of 2 for tubes 1-5 and a dilution  
305 factor of 1 for tube 6. To calculate the true haematocrit values, complete the final *corrected haematocrit for*  
306 *non-haemolysed blood* column in Table 4 by using the following equation:

307

308 **Corrected haematocrit = observed haematocrit x dilution factor**

309

310 *Estimating haemolysis*

311 10. Following the measurement of haematocrit, estimate the percentage of haemolysis of the red blood cells  
312 in the various solutions. To do this, centrifuge the remaining contents of the 6 plastic centrifuge tubes at  
313 6000 rpm for 2 minutes using the centrifuge tube rotor. Take six clean 1.5 ml plastic Eppendorf tubes also  
314 labelled 1-6, pipette 1.5 ml of the supernatant from each correspondingly labelled centrifuge tube into the  
315 labelled Eppendorf tube taking care not to disturb the red blood cell pellet at the bottom of the tube.

316 11. The colours of the six samples of supernatant can then be compared to that of the known haemoglobin  
317 standard solutions prepared at the beginning of the practical. Using the colours of the known haemoglobin  
318 standard solutions as a scale, estimate the concentration of observed supernatant haemoglobin with the  
319 darker the colour of the sample indicating the greater the amount liberated haemoglobin in the supernatant  
320 and hence the greater degree of haemolysis. Using the known haemoglobin concentrations (g/dl) calculated

321 in Table 3, record these observations for tubes 1-6 in the *observed supernatant [Hb]* column of Table 5. If the  
322 colours are not exact matches, estimate whereabouts between the two standards the concentration falls.  
323 12. To convert the observed haemoglobin concentration into an estimated percentage of haemolysis of the  
324 red blood cells, with the exception of the non-haemolysed blood sample which contained 7.5 g/dl  
325 haemoglobin, the blood in the other mixtures was diluted 50:50 and therefore contained half the original  
326 haemoglobin. To estimate the amount of haemolysis that occurred in each sample use the following  
327 calculations and complete the *estimated % haemolysis* column of Table 5:

328

329 **Diluted samples estimated % haemolysis = [observed supernatant Hb concentration / 3.75] x 100**

330 **Non-haemolysed blood estimated % haemolysis = [observed supernatant Hb concentration / 7.5] x 100**

331

332 13. The corrected haematocrit recorded in Table 4 was generated by non-haemolysed red blood cells only as  
333 these were the whole cells that would have made up the packed cell volume in the haematocrit tubes. To  
334 correct for haemolysis in each sample and allow an estimate of what haematocrit would be had there had  
335 been no cell lysis, use the following calculation and complete the final column (*corrected haematocrit*) of  
336 Table 5:

337

338 **Corrected haematocrit (%) = [100 / (100 – estimated % haemolysis)] x corrected haematocrit for non-**  
339 **haemolysed blood (%)**

340

341 **Table 5: Estimated haemolysis and final corrected haematocrit**

342 14. When all data have been collected, each group should pool their final corrected haematocrit (%) data  
343 from Table 5 with the rest of the class using a spreadsheet on a central computer to ensure that group data  
344 can be distributed for more comprehensive analysis following the class.

345

346

347 **Troubleshooting**

348 A common student mistake in this practical is the incorrect or lack of labelling of tubes and pipettes  
349 containing the different solutions during the various steps undertaken. As a result, students lose track of the  
350 contents of tubes they are testing and find their results are meaningless. This is an important error to impress  
351 upon the students as, if this kind of mistake occurs in a clinical setting, the outcome could be life-  
352 threatening. Trained demonstrators should be on hand to spot mistakes early and help students rectify them  
353 as soon as possible.

354

355 Students also often find it difficult to perform the correct calculations to work out dilutions (15). It is  
356 recommended that students are encouraged to attempt these calculations (Tables 1-3) prior to the practical  
357 and come prepared to have these calculations checked by a demonstrator in the practical before proceeding.

358

#### 359 **Safety Considerations**

360 Despite the risk to humans from animal blood being extremely low, when dealing with blood, standard safety  
361 precautions must be taken to minimise the risk of infection. At all times in the laboratory general laboratory  
362 safety rules must be followed including wearing a laboratory coat and using disposable gloves and hand  
363 washing before leaving the laboratory. Any spilled blood or fluids must be wiped up immediately and  
364 disposed of in waste bags provided. All sharps should be disposed of in a sharps box.

365

366 Unless the students are already trained and experienced in using centrifuges, the centrifuge should only be  
367 operated by trained personnel and students should not be left to spin their samples unsupervised. The  
368 centrifuge should be inspected for damage regularly. When using the centrifuge, ensure the tubes are  
369 undamaged, firmly sealed and haven't been overfilled. When placing the tubes in the rotor they must be  
370 balanced and the lid must never be opened while the rotor is moving. The centrifuge should not be left  
371 unattended during use.

372

373

#### 374 **RESULTS**

## 375 **Expected Results**

### 376 *Non-haemolysed blood*

377 Tube 6 which contains the non-haemolysed blood sample prepared in step 5 should be used as a control and  
378 reference point against which to compare any changes to haematocrit in the other blood samples (tubes 1-  
379 5) that were exposed to permeant and non-permeant solutes. Completed sample data tables (Tables 6 and  
380 7) are given here from experiments carried out using pig blood but caution should be taken when making  
381 direct comparisons to the values obtained as, although the relative changes should be the same, the actual  
382 values can vary greatly depending on the blood sample used.

383

### 384 *The effects of hypertonic NaCl*

385 In step 5, non-haemolysed blood was exposed to 2.7% w/v NaCl solution which has an osmolarity of 859  
386 mosmol/l and is hypertonic relative to plasma (tube 1). When red blood cells are placed in a hypertonic  
387 solution, the higher effective osmotic pressure of the bathing solution compared to the intracellular fluid  
388 results in water moving down its osmotic gradient and a net movement of water out of the cell via osmosis  
389 (10). The red blood cells therefore lose their normal biconcave shape and shrink or crenate. This collapse  
390 leads to a decrease in the packed cell volume, or haematocrit, of the solution in comparison to that of the  
391 non-haemolysed blood as the cells take up less space due to the rapid loss of water. Very little haemolysis of  
392 the red blood cells in the solution should be observed as no cells have taken on an additional water load and  
393 burst or haemolysed; however, a few cells may have been damaged during handling and release some  
394 haemoglobin.

395

### 396 *The effects of isotonic NaCl*

397 In step 5, non-haemolysed blood was exposed to an isotonic solution of 0.9% w/v NaCl (osmolarity 286  
398 mosmol/l) (tube 2). This environment has an even distribution of osmolyte particles across both sides of the  
399 cell membrane as intracellular fluid also has an osmolarity around 286 mosmol/l. There is therefore no net  
400 water movement between the bathed red blood cells and the NaCl solution. The haematocrit of the solution



401 should be unaffected and the value similar to that of the non-haemolysed blood. Similarly, little if any  
402 haemolysis of the red blood cells should have occurred.

403

#### 404 *The effects of hypotonic NaCl*

405 In step 5, non-haemolysed blood was exposed to a low osmolarity (143 mosmol/l) hypotonic solution (0.45%  
406 w/v NaCl) (tube 3). When red blood cells are exposed to these conditions where there is a higher  
407 concentration of water and lower effective osmotic pressure outside the cell compared to the intracellular  
408 fluid, this results in net movement of water into the cells via osmosis (11). The cells will increase in size and  
409 some may haemolyse. In this sample therefore a small proportion of haemolysis should have been observed  
410 with increased haemoglobin in the supernatant when compared to the whole blood and the remaining cells  
411 which hadn't lysed would increase in size causing the haematocrit to increase.

412

#### 413 *The effects of distilled water*

414 In step 5, the cells in tube 4 that were bathed in distilled water underwent complete haemolysis and the  
415 estimated % haemolysis should have been 100%. With no ions present in the bathing solution this solution  
416 was very hypotonic resulting in net movement of water into the red blood cells via osmosis causing all the  
417 cells to lose the integrity of their membranes and to haemolyse releasing haemoglobin into the supernatant,  
418 hence the strong red colour of the sample. The resultant corrected haematocrit was 0% as there were no  
419 remaining complete red blood cells to contribute to pack cell volume. Comparing the results of distilled water  
420 (tube 4) and 0.45% w/v (tube 3) is a clear example of how the osmotic fragility or susceptibility of red blood  
421 cells to haemolysis depends on the degree of hypotonicity of the bathing solution.

422

#### 423 *The effects of isosmotic urea*

424 In contrast to NaCl, the membrane is permeable to urea. In Step 5 when red blood cells were bathed in  
425 isosmotic urea (286 mosmol/l) (tube 5), the effects of the permeability of the membrane to urea on both  
426 haematocrit and degree of haemolysis were very different than when red blood cells are exposed to  
427 isosmotic NaCl (tube 2). In the presence of an isosmotic urea solution, the red blood cells underwent

complete haemolysis with a corrected haematocrit of 0%. This is because although isosmotic, the urea solution is not isotonic as urea can freely diffuse across the cell membrane into the cell via passive diffusion and through urea transporters (20, 21). This leads to a change in cell volume as a result of osmotic water movement (13). The isosmotic urea solution is therefore hypotonic because the reflection coefficient of the membrane (permeability) for urea is 0.024 compared to a reflection coefficient of the membrane of 0.3 for NaCl. If the membrane is completely impermeable to a solute the reflection coefficient would be 1. The consequence of this is that the effective osmotic pressure of a urea solution is lower than that of NaCl of the same osmolarity and, as a result, the osmotic gradient across the cell membrane is increased and water moves into the red blood cells via osmosis causing the cell membrane to rupture and the cell to haemolyse. Conversely, NaCl dissociates into  $\text{Na}^+$  and  $\text{Cl}^-$  particles that cannot cross the cell membrane and therefore generate an equal effective osmotic pressure between the extracellular fluid and the intracellular fluid. Under these conditions the osmotic gradient across the cell membrane is maintained and the solution is both isosmotic and isotonic. The same strong red colour of the urea sample in tube 5 should have been observed as that of the distilled water sample in tube 4 as there is 100% haemolysis and 0% corrected haematocrit.

442

#### 443 *Conclusions*

444 The observations and conclusions that should have been drawn from this practical are fundamental to  
445 understanding basic cell physiology. A good grasp of the concepts covered by this practical will help students  
446 appreciate the fact that cell membranes are indeed selectively-permeable and that the tonicity and  
447 osmolarity of fluids affect cell size and structure. This is essential in understanding the concept of  
448 homeostasis and will be referred to in later parts of many physiology courses including during study of the  
449 gastrointestinal tract, regulation of NaCl by the nephron in the renal system and in particular the effect of  
450 dehydration on the whole body.

451

452 Caution must be taken with the practical to ensure students observe the expected results, common mistakes  
453 such as poor labelling of samples and contamination with urea due to pipette confusion can lead to students

454 obtaining results that may not be as anticipated. Careful supervision of students and pooling of data to  
455 analyse the class averages should help prevent this.

456

457 **Table 6: Completed sample results of haematocrit measurements for non-haemolysed blood**

458

459 **Table 7: Completed sample results of estimated haemolysis and final corrected haematocrit**

460

461 **Figure 2: Mean corrected haematocrit values of blood following exposure to permeable and non-**  
462 **permeable solutions (n=64). Error bars represent SEM.**

463

#### 464 **Misconceptions**

465 From early on in many physiology-based courses, students struggle with the concepts of osmolarity and  
466 tonicity and find it difficult to relate them to the direction of water movement. This practical can help  
467 students to visualise different solutions and the effects that these can have on red blood cells. By being able  
468 to see the colour change of the non-haemolysed blood samples mixed with various solutions they can relate  
469 the theory of osmosis to what has happened to the cells when water moves out of the red blood cells in a  
470 hypertonic solution (tube 1), when there is no net movement of water in an isotonic solution (tube 2) and  
471 when water moves into the cells in a hypotonic solution (tube 3). The practical also helps students with  
472 misconceptions surrounding haemolysis and haematocrit. The students often mistake the contents of the  
473 tubes (4 and 5) containing distilled water and urea as having 100% haematocrit due to the dark colouration  
474 of the whole sample with little or no visible plasma band rather than the product of 100% haemolysis due to  
475 the effects of the tonicity. Trained demonstrators should be on hand to ensure that students are able to  
476 relate their findings to the learning outcomes and a whole class tutorial on the outcomes of the experiments  
477 and their meanings is scheduled for one week after the practical class.

478

#### 479 **Evaluation of Student Work**

480 As discussed, there are a number of ways in which students may not get the results expected. To ensure that  
481 students have access to some representative data on which to perform any post-practical analysis, students  
482 are expected to pool their results with the rest of the group to produce group data for the class prior to the  
483 finish and this is then shared online with the students to use immediately after the class.

484

485 We assess our students on this work through the submission of an online post-practical assessment. This  
486 takes the form of multiple choice questions on the background physiology, method, results and physiological  
487 significance of the findings. Questions included in the assessment test whether students have grasped the  
488 direction of fluid movement in the presence of different solutions and the physiological reasons for this. With  
489 regards to further data analysis, students are expected to plot the class findings in graphical form. They are  
490 instructed to produce a column/bar graph of the final corrected haematocrits (%) from all solutions with the  
491 error bars as the standard error of the mean and this is uploaded as part of their online post-practical  
492 assessment (see Figure 2). Students are also given information on how to calculate osmolarity and are  
493 expected to perform calculations themselves.

494

#### 495 **Inquiry Applications**

496 As this practical is run at the start of the students' exploration of physiology at undergraduate level there is  
497 very limited inquiry in this practical and would be considered 'Methods' level. The questions being explored  
498 and the procedure being followed are clearly set out by the teacher running the practical. Students carry out  
499 the practical, with assistance from demonstrators and analyse the data during the session.

500

501 However, there is scope for this practical to become more student-centred and to be used at a higher inquiry  
502 level by making a number of modifications to the protocol. These could include giving the students less rigid  
503 instructions on how and what haemoglobin standard solutions to produce and allowing them to make a wider  
504 range of standard solutions to allow more accuracy in estimating the final degree of haemolysis. They could  
505 also be asked to predict the impact of the different solutions on the haematocrit and haemoglobin  
506 concentration and subsequently test these predictions. More sophisticated techniques could also be used to

507 measure haemolysis in the supernatant such as spectrophotometry and using light microscopy to visualise  
508 the red blood cells after they have undergone crenation or haemolysis as opposed to the estimates made in  
509 this experiment by eye.

510

511 This practical could be incorporated into a number of different biomedical programmes, from  
512 Biology/Physiology/Biochemistry honours programmes to aid in the understanding of the fundamental  
513 concepts of cell transport and membrane structure as well as developing vital scientific skills including  
514 handling blood and performing serial dilutions. Professional disciplines such as medicine and veterinary  
515 science would also benefit from this practical to further explore the concepts of osmotic fragility and the  
516 administration of intravenous fluids and the clinical implications that these can have.

517

518 The concept of osmotic fragility could be explored further by using a series of hypotonic solutions and  
519 recording % haemolysis. From these data an osmotic fragility curve could be plotted to explore the internal  
520 pressures exerted on the cell membrane when water diffuses into a cell. Students could further consider how  
521 the shape of red blood cells e.g. sickle cells may affect the osmotic fragility of red blood cells (8) and how  
522 haemolytic diseases such as thalassemia and hereditary spherocytosis are a result of changes in osmotic  
523 fragility both by extending this practical and the use of further resources (12).

524

525 Further examples of the differences between tonicity and osmolarity and the effect of permeable solutes can  
526 be used, for instance with the addition of glucose, a particularly clinically relevant solute with regards to  
527 intravenous administration of fluids. The clinical application of an understanding of this concept can be  
528 emphasised including those surrounding patient safety (17).

529

#### 530 **ADDITIONAL RESOURCES**

531 For additional information on this topic, any undergraduate level physiology textbook should provide  
532 relevant background information required to understand the theory this practical is based upon.

533

534   **ACKNOWLEDGEMENTS**

535   We would like to thank Mr. Dave Gee, Senior Physiology Teaching Laboratory Technician and his technical  
536   team for the invaluable help they provide in preparing this practical and for their assistance with this article.

537   We acknowledge the contribution of Professor Graham Mitchell for the initial design of this experiment.

538

539   **REFERENCES**

- 540       1. Bird GW. The red cell. *British Medical Journal*. 1(5795):293-7, 1972.
- 541       2. Choong K, Kho ME, Menon K, Bohn D. Hypotonic versus isotonic saline in hospitalised children: a  
542       systematic review. *Archives of Disease in Childhood*. 91(10):828-35, 2006.
- 543       3. Cooper GM. *The Cell: A Molecular Approach*. 6<sup>th</sup> edition. Sunderland (MA): Sinauer Associates, 2013.
- 544       4. Dourmashkin RR, Rosse WF. Morphologic changes in the membranes of red blood cells  
545       undergoing hemolysis. *American Journal of Medicine*. 41(5):699-710, 1966.
- 546       5. Feher JJ, Ford GD. A simple student laboratory on osmotic flow, osmotic pressure, and the reflection  
547       coefficient. *Advances in Physiology Education*. 13(1): S10-S20, 1995.
- 548       6. Gorter E, Grendel F. On bimolecular layers of lipids on the chromocytes of the blood. *Journal of*  
549       *Experimental Medicine*. 41:439-443, 1925.
- 550       7. Hasler U, Leroy V, Martin PY, Féraille E. Aquaporin-2 abundance in the renal collecting duct: new  
551       insights from cultured cell models. *American Journal of Physiology Renal Physiology*. 297(1):F10-8,  
552       2009.
- 553       8. Hijiya N, Horiuchi K, Asakura T. Morphology of sickle cells produced in solutions of varying  
554       osmolarities. *Journal of Laboratory and Clinical Medicine*. 117(1):60-6, 1991.
- 555       9. Khan Academy, Structure of the Plasma Membrane.  
556       [www.khanacademy.org/science/biology/membranes-and-transport/the-plasma-](http://www.khanacademy.org/science/biology/membranes-and-transport/the-plasma-membrane/a/structure-of-the-plasma-membrane)  
557       [membrane/a/structure-of-the-plasma-membrane](http://www.khanacademy.org/science/biology/membranes-and-transport/the-plasma-membrane/a/structure-of-the-plasma-membrane) Accessed 04/11/16.
- 558       10. Kregenow FM. The response of duck erythrocytes to hypertonic media. Further evidence for a  
559       volume-controlling mechanism. *Journal of General Physiology*. 58(4):396-412, 1971a.

11. Kregenow FM. The response of duck erythrocytes to nonhemolytic hypotonic media. Evidence for a volume-controlling mechanism. *Journal of General Physiology*. 58(4):372-95, 1971b.
12. Kumar S. An analogy for explaining erythrocyte fragility: concepts made easy. *Advances in Physiology Education*. 26(2): 134-135, 2002.
13. Lezama R, Díaz-Téllez A, Ramos-Mandujano G, Oropeza L, Pasantes-Morales H. Epidermal growth factor receptor is a common element in the signaling pathways activated by cell volume changes in isosmotic, hyposmotic or hyperosmotic conditions. *Neurochemical Research*. 30(12):1589-97, 2005.
14. Macey RI. Transport of water and urea in red blood cells. *American Journal of Physiology*. 246(3 Pt 1):C195-203, 1984.
15. Moni RW, Hryciw DH, Poronnik P, Lluca LJ, Moni KB. Assessing core manipulative skills in a large, first-year laboratory. *Advances in Physiology Education*. 31(3): 266-269, 2007.
16. Preston GM, Carroll TP, Guggino WB, Agre P. Appearance of water channels in *Xenopus* oocytes expressing red cell CHIP28 protein. *Science*. 256(5055): 385-7, 1992.
17. Silverthorn D. 2010. Using IV fluid therapy to teach the principles of osmolarity & tonicity. MedEdPORTAL, <http://services.aamc.org/30/mededportal/servlet/s/segment/mededportal/?subid=8082>. Accessed 06/05/16.
18. Singer SJ, Nicolson GL. The Fluid Mosaic Model of the Structure of Cell Membranes. *Science*. 175(4023): 720-731, 1972.
19. Strange K. Cellular volume homeostasis. *Advances in Physiology Education*. 28(4): 155-159, 2004.
20. Wieth JO, Funder J, Gunn RB, Brahm J. Passive transport pathways for chloride and urea through the red cell membrane. In: *Comparative Biochemistry and Physiology of Transport*, edited by Bolis K, Bloch K, Luria SE, and Lynen F. Amsterdam: Elsevier/North-Holland. 1974, 317-337.
21. Yang B. Transport characteristics of urea transporter-B. *Subcellular Biochemistry*. 73:127-35, 2014.